



QUALITY AT IT'S BEST..

MANIFOLD






ARKA - THE PERFECT MANIFOLD

Increased usability with smoothly rotating handles is particular to meet all procedural needs



High Pressure Stopcock



CATALOG NO	Product
A1113	500 PSI High Pressure Stopcock

CATALOG NO	Product
A1069	Manifold Off Right
A3841	Manifold Off Left
A1076	Manifold Set
A1014	Angiography Syringe 12 ml

Pressure: 500 PSI, holding higher pressure than most competitors
Innovative ergonomic handles - smooth turning torque design improves ease of use and convenience.

Elevated handle design enables turning of the handles without touching the ports

High profile design for easy cleaning and disinfecting

Transparent housing for visual control

Available with 2 or 3 ports, OFF or ON handles, Right or Left orientation, different handle colors including color coded handles,

