



QUALITY AT IT'S BEST..

RADIAL INTRODUCER SHEATH






ARKA - RADIAL INTRODUCER SHEATH

FEATURES & BENEFITS

- Thin **radiopaque sheath** with **anti-kinking sleeve** for easy catheter handling
- Dilator is **tapered** - **minimizes vessel trauma**
- **Wire type** - Stainless steel wire, plastic Jacketed with Hydrophilic coating
- Mindfully designed **Haemostasis Valve** to avoid blood reflux , air aspiration
- **Suture ring** on hub convenient for securing the sheath in multi-stage procedures
- **High quality raw materials** and expertise in precision insertion offers a smooth sheath and dilator for an simplified
- Available as **Screw lock** and **push-lock** designs.



Mindfully designed valve assures reliable homeostasis and facilitates catheter manipulation remarkably.



Our precision expertise improves dilator and sheath transitions, ensuring smooth connections and a dual tapered design that reduces trauma risk to surrounding tissue at the puncture site.



3-way stopcock at the tip of extension tube allows for proper flush of dilator and assures clear visualization of bleeding

ORDERING INFORMATION

PRODUCT CODE	NEEDLE TYPE	NEEDLE LENGTH (CM)	NEEDLE GAUGE	FR SIZE	INTRODUCER LENGTH (CM)	INTRODUCER TOP DIAMETER (MM)	TIP TYPE	WIRE TYPE	WIRE LENGTH (CM)	WIRE DIAMETER (MM)	BOX QTY
IS9011	Radial Arterial	4	21 ga	4	7	0.021	Straight	Stainless steel wire, plastic Jacketed with Hydrophilic coating	80	0.021	5
IS9012				5	7	0.021			50		
IS9013				6	7	0.021			50		
IS9016				5	11	0.021			80		
IS9017				6	11	0.021			80		

