



QUALITY AT IT'S BEST..


# ARKA RADIAL ANGIOGRAPHY CATHETER





# ARKA - RADIAL ANGIOGRAPHY CATHETER

The ARKA Diagnostic Catheter Portfolio features high-density braiding for excellent responsiveness, offering dedicated and universal nylon catheters with superior pushability and kink resistance. The SUPER TORQUE PLUS portfolio, made of polyurethane, provides one-to-one torque control, kink resistance, and consistent lumen diameter across various shapes.



Soft tip reduces potential for vessel trauma  
Double braided stainless steel for 1:1 Torque control  
Maximum Guide wire: 0.038" (0.97 mm.)  
Maximum Pressure: 1200psi (84 bar)

## ORDERING INFORMATION

AC-TIG	CATETER CARDIOLOGICO PERF ULTIMATE 1 RADIAL 5FR 100CM 038	ACTIG15F
	CATETER CARDIOLOGICO PERF ULTIMATE 2 RADIAL 5FR 100CM 038	ACTIG25F
	CATETER CARDIOLOGICO PERF ULTIMATE 1 RADIAL 6FR 100CM 038	ACTIG16F