



# Arrivaa<sup>TM</sup> SC

## Balloon Dilatation Catheter



**Advanced Laser welding Technology** to realize seamless connection between head and tube.

**Tri-Fold Shape Memory Technology** offers excellent rewarpability for multiple inflation uses.

**Silicone Resin into The Guidewire Lumen** to create a smooth, low-friction path for seamless guidewire movement through the catheter.



# Arrivaa<sup>TM</sup> SC

## Balloon Dilatation Catheter



### TECHNICAL SPECIFICATIONS

|                         |                         |
|-------------------------|-------------------------|
| Proximal Shaft          | 1.90 F                  |
| Distal Shaft            | 2.60 F                  |
| Catheter Working Length | 140 cm                  |
| Tip Length              | 3.0 mm tip for 01.0-4.0 |
| Marker Bands            | 2                       |
| Coating                 | BSL-Hydrophilic         |
| Nominal Pressure        | 6 ATM                   |
| Rated Burst Pressure    | 14 ATM                  |
| Folds                   | 3                       |

### ORDERING INFORMATION

| DIA<br>↓ | Length<br>↗ | LENGTH |    |    |    |
|----------|-------------|--------|----|----|----|
| 0.75     | 05          | 10     |    |    |    |
| 1.00     | 08          | 10     | 12 | 15 | 20 |
| 1.25     | 08          | 10     | 12 | 15 | 20 |
| 1.50     | 08          | 10     | 12 | 15 | 20 |
| 1.75     | 08          | 10     | 12 | 15 | 20 |
| 2.00     | 08          | 10     | 12 | 15 | 20 |
| 2.50     | 08          | 10     | 12 | 15 | 20 |



Marketed by:  
**Zon Medical Devices**  
QUALITY AT ITS BEST.

The Netherland Office  
4-15 EdelBox Osloweg  
4 9723 BL Groningen  
The Netherlands

Hyderabad Office

9 No. 2-3-125, Plot No. 36/B,  
Second Floor, Road No. 8,  
Co-operative Bank Colony,  
Nagole, Hyderabad T.S.,  
India 500068.

Noida Office

9 5th Floor, B 522, Logix  
Technova Block B,  
Goutambudh Nagar, Noida  
Expressway, Noida Sector  
132 - 201304