




# FEMORAL ANGIOGRAPHY CATHETER

AVAILABLE AS OEM/ODM



The Diagnostic Catheter Portfolio features high-density braiding for excellent responsiveness, offering dedicated and universal nylon catheters with superior pushability and kink resistance. The SUPER TORQUE PLUS portfolio, made of polyurethane, provides one-to-one torque control, kink resistance, and consistent lumen diameter across various shapes.



Soft tip reduces potential for vessel trauma

Double braided stainless steel for 1:1 Torque control

Maximum Guide wire: 0.038" (0.97 mm.)

Maximum Pressure: 1200psi (84 bar)

# ORDERING INFORMATION

	TIP SHAPE	Spesifikasyon	CATALOG NO
ANGIOGRAPHY CATHETER	AC-JL	CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JL3.5	ACJL354F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL3.5	ACJL355F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL3.5	ACJL356F
		CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JL4.0	ACJL404F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL4.0	ACJL405F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL4.0	ACJL406F
		CATETER CARD JUDKINS HIGHFLOW 5FR X 100CM 038 JL4.0	HACJL405F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL4.5	ACJL455F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL4.5	ACJL456F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL5.0	ACJL505F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL5.0	ACJL506F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL6.0	ACJL605F
	AC-JR	CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JR3.5	ACJR354F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR3.5	ACJR355F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR3.5	ACJR356F
		CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JR4.0	ACJR404F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR4.0	ACJR405F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR4.0	ACJR406F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR5.0	ACJR505F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR5.0	ACJR506F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR6.0	ACJR605F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR6.0	ACJR606F
	AC-AL	CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL1	ACAL15F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL1	ACAL16F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL2	ACAL25F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL2	ACAL26F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL3	ACAL35F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL3	ACAL36F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL4	ACAL45F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL4	ACAL46F
	AC-AR	CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR1	ACAR15F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR1	ACAR16F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR2	ACAR25F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR2	ACAR26F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR3	ACAR35F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR3	ACAR36F
	AC-PIG	CATETER PIGTAIL ULTRA HIGH FLOW 5FR X 65CM	ACPIGS5F1065
		CATETER CARD PIGTAIL 145° 6FR X 110CM 038	ACPIGA6F
		CATETER CARD PIGTAIL 145° 5FR X 110CM 038	ACPIGASF
		CATETER CARD PIGTAIL RETO 4FR X 80CM 038	ACPIGS4F1080
		CATETER CARD PIGTAIL RETO 4FR X 110CM 038	ACPIGS4F
		CATETER CARD PIGTAIL RETO 5FR X 110CM 038	ACPIG5F
		CATETER CARD PIGTAIL RETO 5FR X 100CM 038 20 MARKER BANDS	ACPIG5F20M10H
		CATETER CARD PIGTAIL RETO 5FR X 125CM 038	ACPIG5F
		CATETER CARD PIGTAIL RETO 6FR X 110CM 038	ACPIG6F
		CATETER PIT TALL 5FR X 100CM 038 20BANDS(1CM APART) 5 SH	ACPIGT2015F
		CATETER PIT TALL 5FR X 100CM 038 20BANDS(1CM APART) 10 SH	ACPIGT2025F
		CATETER PIT TALL 5FR X 65CM 038 11BANDS(1CM APART) 10 SH	ACPIGT1115F165
		CATETER PIT TALL 5FR X 100CM 038 11BANDS(1CM APART) 10 SH	ACPIGT1115F
		CATETER PIT TALL 5FR X 65CM 038 2 MARKER BANDS(2CM APART) 10 SH	ACPIGT215F165
	AC-LCB	CATETER CORONARIO BYPASS LEFT 5FR X 100CM 038 LCB	ACLCB5F
		CATETER CORONARIO BYPASS RIGHT 6FR X 100CM 038 RCB	ACRCB6F
		CATETER CORONARIO BYPASS RIGHT 5FR X 100CM 038 RCB	ACRCB5F
	AC-IM	CATETER MAMARIA INTERNA PERFORMA 5FR X 100 CM	ACIM5F
	AC-CIM	CATETER MAMARIA INTERNA PERFORMA 4FR X 100 CM 038	ACIM4F
		CATETER MAMARIA INTERNA PERFORMA 6FR X 100 CM	ACIM6F
	AC-MPA	CATETER CARD MULTIPURPOSE 4FR X 100CM 038 A2	ACMPA24F
		CATETER CARD MULTIPURPOSE 5FR X 100CM 038 B2	ACMPA25F
		CATETER CARD MULTIPURPOSE 4FR X 100CM 038 A1	ACMPA14F
		CATETER CARD MULTIPURPOSE 6FR X 100CM 038 A1	ACMPA16F
		CATETER CARD MULTIPURPOSE 5FR X 100CM 038 A1	ACMPA15F
		CATETER CARD MULTIPURPOSE 6FR X 100CM 038 A2	ACMPA26F