




FEMORAL INTRODUCER SHEATH

AVAILABLE AS OEM/ODM






- Thin radiopaque sheath with anti-kinking sleeve for easy catheter handling
- Dilator is tapered - minimizes vessel trauma
- Wire type - Stainless steel wire, plastic Jacketed with Hydrophilic coating
- Mindfully designed Haemostasis Valve to avoid blood reflux , air aspiration
- Suture ring on hub convenient for securing the sheath in multi-stage procedures
- High quality raw materials and expertise in precision insertion offers a smooth sheath and dilator for an simplified
- Available as Screw lock and push-lock designs.



FEATURES & BENEFITS



Mindfully designed valve assures reliable haemostasis and facilitates catheter manipulation rem bly.



Our expertise in precision brings the dilator and sheath transition to an excellent level. The smooth between components and dual tapered design minimizes the risk of trauma in surrounding tissue of puncture point.



3-way stopcock at the tip of extension tube allows for proper flush of dilator and assures clear visualization of bleeding

SIZES

Available in 4 5 6 7 8 9 10

ACTUAL SIZE	FRENCH	MM	INCHES
●	4	1.35	0.053
●	5	1.67	0.066
●	6	2	0.079
●	7	2.3	0.092
●	8	2.7	0.105
●	9	3	0.118
●	10	3.3	0.131



ORDERING INFORMATION

PRODUCT CODE	NEEDLE TYPE	FR SIZE	INTRODUCER LENGTH (CM)	INTRODUCER TOP DIAMETER (MM)	TIP TYPE	WIRE TYPE	WIRE LENGTH (CM)	WIRE DIAMETER (MM)	BOX QTY
IS9000	18G 7CM	4	9	0.035	J-curved	Stainless steel wire, plastic Jacketed with Hydrophilic coating	50	0.035	10
IS9001		5	11	0.038			45	0.038	10
IS9002		6	11	0.038			45	0.038	10
IS9003		7	11	0.038			45	0.038	10
IS9004		8	11	0.038			45	0.038	10
IS9005		9	11	0.038			45	0.038	10
IS9006		10	11	0.038			45	0.038	10
IS9007		6	23	0.038			80	0.038	10
IS9008		7	23	0.038			80	0.038	10
IS9009		8	23	0.038			80	0.038	10

