





HYDROPHILIC GUIDE WIRE

AVAILABLE AS OEM/ODM




HYDROPHILIC GUIDEWIRE FOR PERFECTION



WHY GUIDEWIRES?

- Significance** In every coronary and peripheral procedure there comes a significant action to take: gain access to the treatment zone. This is why guidewires are essential actors in vascular interventions and this is why offers "the best".
- Core** guidewires are made of excellent quality torque-stable cores with optimal maneuverability and enhanced kink resistance particular for effective and atraumatic navigation through challenging vascular anatomies.
- Coating** Available with Hydrophilic coatings, guidewires feature a smooth surface to ensure rapid movement, secure advancement and device placement.
- Tip Configuration** guidewires are available with straight, angled and j-curved tip configurations. The proximal shaft with the ideal stiffness assures perfect pushability while the distal tip enables the operator for smooth and precise advancement.



TECHNICAL SPECIFICATIONS

Usable length	150 cm, 180 cm, 260 cm
Diameter	0,035" (0.89mm)
Core Material	Nitinol for Hydrophilic coating
Coating	Hydrophilic
Sterilization	Sterilized with ethylene oxide gas
Tip Configuration	Straight, Angled
Packaging	Dispenser and Tyvek package

ORDER INFORMATION

Ref code	Coating	Tip	Diameter	Length
GW2205	Hdyrophilic	Angled/Straight	0,035	150 cm
GW2206	Hdyrophilic	Angled/Straight	0,035	180 cm
GW2207	Hdyrophilic	Angled/Straight	0,035	260 cm