



INFLATION DEVICES

AVAILABLE AS OEM/ODM



Inflation Device

Effortless Inflation
Successful Operation



Sidra

One-Hand One-Touch

Effort-Free Error-Free

- One-hand use for easy and quick inflation and deflation
- The rotating handle and ergonomic device body is particular for comfortable grip and handling
- Features a 30 ATM pressure rating with a 30 mL syringe volume.
- The pressure gauge is available with 35 ATM option.



Classic

Perceivable Power



Gun Model

Remarkable Design

CATALOG NO	Product
A1052	Inflation Device 20ml, 30atm Classic Type
A2967	Inflation Device 30ml, 30atm Classic Type
A3780	Inflation Device 20ml, 30atm Gun Type
A3797	Inflation Device 30ml, 30atm Ergonomic Type
A2967	Ptca Kit (Balloon Inflation Device + Hemostasis Valve Set)